

Quattro™ Air ...
the lighter way
to better sleep.



ResMed's new lightweight **Quattro Air** Full Face Mask builds on Quattro's proven performance and reliability. It provides enhanced comfort and ease of use, and is 45% lighter* than the previous Mirage Quattro mask.

Benefits for you ...

Builds on the strength of Quattro

Quattro Air provides performance you can rely on, in a new sleeker, lighter mask with less bulk on the face, so it's bound to have strong patient appeal.

Less parts for quick and easy set-ups

Quattro Air has fewer parts than most other Full Face Masks (FFM) – just four. It's quick to demonstrate and easy to train patients to use, supporting patient compliance.

Improved patient usability

The mask features greatly reduced complexity with improved comfort and stability, reducing patient call-backs and freeing up more of your time.

Innovative design for easy acceptance

It's light, streamlined and comfortable, allowing you to assist those of your patients who struggle to stick with therapy.

... and your patients

Lightweight and less bulk on the face

Quattro Air is light with a clearer field of vision and less bulk on the face than other FFM's, so patients are less likely to feel claustrophobic or intimidated.

Intuitive and uncomplicated

Quattro Air is simple for patients to use, for quick and easy cleaning and reassembly.

Forehead stability

Its advanced Flex-wing forehead support is comfortable and feels light on the face.

Quiet operation

The circular venting gently and quietly diffuses airflow away for less user and partner disturbance—only 28 dBA.**

* The Medium size of Quattro Air is 45% lighter than the Medium size of Mirage Quattro.

** Declared dual numbered noise emission values (with uncertainty 3 dBA) in accordance with ISO 4871.e1.



Quattro™ Air
FULL FACE MASK

Even more comfort ...

Quattro Air is lighter, sleeker, and even more comfortable than the previous Mirage Quattro mask due to a series of design enhancements.

It's also simpler to use, trimming down the number of parts to just four. That makes it easier for patients to understand how to use their Quattro Air and simple for them to put it back together after cleaning.

SoftEdge™ headgear design

has premium rolled-edge fabric that reduces facial marking and is also comfortable to wear



Improved comfort

and less chance for skin irritation at the sensitive nasal bridge area

Circular diffused venting

directs air away for less user and partner disturbance

Quick release elbow

with squeeze-tabs to make attaching and detaching quicker and easier

Flex-wing forehead support


allows a dynamic range of fit and adjustment, and feels light on the face

Dual wall Spring Air™ cushion

delivers a high-performance seal with even more comfort

Headgear loops

simplify attachment for improved usability



Product Codes

| | |
|---|-------|
| Quattro Air mask system Small | 62701 |
| Quattro Air mask system Medium | 62702 |
| Quattro Air mask system Large | 62703 |
| Quattro Air for Her mask system X-Small | 62740 |
| Quattro Air for Her mask system Small | 62741 |
| Quattro Air for Her mask system Medium | 62742 |

Quattro™ Air
FULL FACE MASK

Quattro™ Air for Her
FULL FACE MASK



ResMed Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA. **EC/REP** ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See www.resmed.com for other ResMed locations worldwide. Quattro, SoftEdge and Spring Air are trademarks of ResMed Ltd and are registered in the U.S. Patent and Trademark Office. © 2013 ResMed Ltd. Specifications may change without notice. 1015994/1 2013-03